

Goat Anti-Mouse IgG(H+L)(Biotin conjugated)

Catalog No. E-AB-1024

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	Biotin
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

IHC	1:100-200
WB	1:300-10000(Enhanced chemiluminescent detection)
WB	1:200-400(DAB detection)
Direct ELISA	1:5000-20000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.