

# Goat Anti-Rabbit IgG(H+L)(Biotin conjugated)

Catalog Number:E-AB-1025 4 Publications



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Rabbit
<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Biotin
<b>Formulation</b>	PBS with 50% glycerol, PH7.4

## Applications Recommended Dilution

<b>Immunohistochemistry</b>	1:100-200
<b>Western Blotting</b>	1:300-10000(Enhanced chemiluminescent detection)
<b>Western Blotting</b>	1:200-400(DAB detection)
<b>Direct ELISA</b>	1:5000-20000(TMB detection)

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017