

Mouse Anti-Rabbit IgG(H+L)(Biotin conjugated)

Catalog No. E-AB-1030

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	Rabbit IgG
Host	Mouse
Purification	Affinity purification
Conjugation	Biotin
Buffer	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Immunohistochemistry	1:100-200
Western Blotting	1:300-10000(Enhanced chemiluminescent detection)
Western Blotting	1:200-400(DAB detection)
Direct ELISA	1:5000-20000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only