

## Goat Anti-Human IgG(H+L)(Biotin conjugated)

**Catalog No.** E-AB-1031

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Reactivity</b>	Human
<b>Immunogen</b>	Human IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Biotin
<b>Buffer</b>	Biotin:Ab=20:1, 0.01M PBS with 50% glycerol, PH7.6

### Applications Recommended Dilution

<b>Immunohistochemistry</b>	1:100-200
<b>Western Blotting</b>	1:300-10000(Enhanced chemiluminescent detection)
<b>Western Blotting</b>	1:200-400(DAB detection)
<b>Direct ELISA</b>	1:5000-20000(TMB detection)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.