

Rabbit Anti-Mouse IgG

Catalog Number: E-AB-1033 1 Publications



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Rabbit
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS, PH7.4

Applications Recommended Dilution

Other	When the antibody is labeled, for example, with HRP, FITC, BIOTIN, etc, BSA and NaN ₃ cannot be contained in reagents.
--------------	---

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017