

Goat Anti-Mouse IgG

Catalog No. E-AB-1035

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Other	When the antibody is labeled, for example, with HRP, FITC, BIOTIN, etc, BSA and NaN ₃ cannot be contained in reagents.
--------------	---

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only