

Rabbit Anti-Goat IgG

Catalog No. E-AB-1038

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Goat
Immunogen	Goat IgG
Host	Rabbit
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS, PH7.4

Applications Recommended Dilution

Other When the antibody is labeled, for example, with HRP, FITC, BIOTIN, etc, BSA and NaN₃ cannot be contained in reagents.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only