

Goat Anti-Rat IgG(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1041

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rat
Immunogen	Rat IgG
Host	Goat
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

WB	1:2000-5000 (Enhanced chemiluminescent detection)
WB	1:300-1500(DAB detection)
Direct ELISA	1:3000-5000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only