

# Streptavidin(FITC conjugated)

Catalog Number:E-AB-1042



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

**Conjugation** FITC (Ex/Em:495nm/525nm)  
**Formulation** 0.01M PBS with 0.05% Proclin300, 1% protective protein, PH7.4

## Applications Recommended Dilution

**FCM** 1:20-1:500  
**IHC** 1:50-1:500  
**ICC** 1:50-1:500

## Preparation & Storage

**Storage** Store at 2~8°C. Avoid repeated freeze / thaw cycles. Protected from prolonged exposure to light.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017