

Streptavidin(peroxidase/HRP conjugated)

Catalog No. E-AB-1043

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

IHC-F	1:100-500
IHC-P	1:100-200
WB	1:3000-10000(Enhanced chemiluminescent detection)
WB	1:500-3000(DAB detection)
Direct ELISA	1:10000-100000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---