

## Avidin(peroxidase/HRP conjugated)

**Catalog No.** E-AB-1044

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

|                     |   |
|---------------------|---|
| <b>Purification</b> | Affinity purification   |
| <b>Conjugation</b>  | HRP   |
| <b>Buffer</b>       | PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4 |

### Applications Recommended Dilution

|   |  |
|---|--|
| <b>Immunohistochemistry in frozen section</b>   | 1:100-500  |
| <b>Immunohistochemistry in paraffin section</b> | 1:100-200  |
| <b>Western Blotting</b>                         | 1:2000-4000(Enhanced chemiluminescent detection) |
| <b>Western Blotting</b>                         | 1:200-400(DAB detection)                         |
| <b>Direct ELISA</b>                             | 1:10000-100000(TMB detection)                    |

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### For Research Use Only