

Goat Anti-Mouse IgG(H+L)(Alkaline Phosphatase/AP conjugated)

Catalog No. E-AB-1045

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	AP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Western Blotting	1:1000-3000(Enhanced chemiluminescent detection)
Direct ELISA	1:1000-3000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---