

## Goat Anti-Rat IgG(H+L)(Alkaline Phosphatase/AP conjugated)

Catalog No. E-AB-1048

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Reactivity</b>	Rat
<b>Immunogen</b>	Rat IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	AP
<b>Buffer</b>	PBS with 50% glycerol, PH7.4

### Applications Recommended Dilution

**Western Blotting:**  
1:1000-3000(Enhanced chemiluminescent detection ) Direct

**ELISA:**  
1:1000-3000(TMB detection)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### For Research Use Only