

Donkey Anti-Goat IgG(H+L)(Alkaline Phosphatase/AP conjugated)

Catalog No. E-AB-1049

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Goat
Immunogen	Goat IgG
Host	Donkey
Purification	Affinity purification
Conjugation	AP
Buffer	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Western Blotting:
1:1000-3000(Enhanced chemiluminescent detection) Direct
ELISA:
1:1000-3000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only