

Donkey Anti-Goat IgG(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1050

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Goat
Immunogen	Goat IgG
Host	Donkey
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Western Blotting	1:2000-5000(Enhanced chemiluminescent detection)
Western Blotting	1:300-1500(DAB detection)
Direct ELISA	1:3000-5000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.