

## Donkey Anti-Mouse IgG(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1051

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Donkey
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Buffer</b>	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

### Applications Recommended Dilution

<b>WB</b>	1:2000-5000 (Enhanced chemiluminescent detection)
<b>WB</b>	1:300-1500(DAB detection)
<b>Direct ELISA</b>	1:3000-5000(TMB detection)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.