

Goat Anti-Mouse IgG(H+L)(Elab Fluor® 594 conjugated)

Catalog No. E-AB-1059

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	Elab Fluor®594 (Ex/Em:591nm/614nm)
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

IF	1:100-300
-----------	-----------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

For Research Use Only