

## Goat Anti-Rabbit IgG(H+L)(Elab Fluor® 594 conjugated)

**Catalog No.** E-AB-1060

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Reactivity</b>	Rabbit
<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Elab Fluor®594 (Ex/Em:591nm/614nm)
<b>Buffer</b>	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

### Applications Recommended Dilution

<b>IF</b>	1:100-300
-----------	-----------

### Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

### For Research Use Only