

Rabbit Anti-Chicken IgY++(IgG)(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1052

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Chicken
Immunogen	Chicken IgY++
Host	Rabbit
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

WB	1:2000-10000(Enhanced chemiluminescent detection)
Direct ELISA	1:1000-20000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---