

Rabbit Anti-Chicken IgY++(IgG)(H+L)(peroxidase/HRP conjugated)



Catalog Number:E-AB-1052

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Chicken
Immunogen	Chicken IgY++
Host	Rabbit
Purification	Affinity purification
Conjugation	HRP
Formulation	PBS with 0.03% Proclin300 and 50% glycerol, PH7.4

Applications Recommended Dilution

Western Blotting	1:2000-10000(Enhanced chemiluminescent detection)
Direct ELISA	1:1000-20000(TMB detection)

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017