

Monoclonal Mouse anti-Human IgE(C3-domain specific)(peroxidase/HRP conjugated)

Catalog No. E-AB-1074

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human
Immunogen	Human IgE C3 domain
Host	Mouse
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Direct ELISA	1:200-500(TMB detection)
---------------------	--------------------------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---