

Goat Anti-Mouse IgG(H+L)(Rhodamine/TRITC conjugated)

Catalog No. E-AB-1064

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	TRITC (Ex/Em:550nm/620nm)
Buffer	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Immunofluorescence:
1:20-100

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.