

Goat Anti-Mouse IgM(H+L)(FITC conjugated)

Catalog No. E-AB-1067

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

| | |
|---------------------|---|
| Reactivity | Mouse |
| Immunogen | Mouse IgM |
| Host | Goat |
| Purification | Affinity purification |
| Conjugation | FITC (Ex/Em:495nm/525nm) |
| Buffer | PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4 |

Applications Recommended Dilution

| | |
|-----------|----------|
| IF | 1:20-100 |
|-----------|----------|

Preparation & Storage

| | |
|----------------|---|
| Storage | Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
|----------------|---|

For Research Use Only