

Goat Anti-Human IgM(H+L)(FITC conjugated)

Catalog No. E-AB-1068

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Human
Immunogen	Human IgM
Host	Goat
Purification	Affinity purification
Conjugation	FITC (Ex/Em:495nm/525nm)
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

IF	1:20-100
-----------	----------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---