

## Goat Anti-Human IgM(H+L)(Biotin conjugated)

Catalog No. E-AB-1069

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Reactivity</b>	Human
<b>Immunogen</b>	Human IgM
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Biotin
<b>Buffer</b>	PBS with 50% glycerol, PH7.4

### Applications Recommended Dilution

**Immunohistochemistry:** 1:100-200  
**Western Blotting:** 1:200-4000(Enhanced chemiluminescent detection )  
**Western Blotting:** 1:200-400(DAB detection)  
**Direct ELISA:** 1:1000-5000(TMB detection)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### For Research Use Only