

Goat Anti-Mouse IgM(H+L)(Biotin conjugated)

Catalog No. E-AB-1070

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgM
Host	Goat
Purification	Affinity purification
Conjugation	Biotin
Buffer	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Immunohistochemistry: 1:100-200
Western Blotting: 1:200-4000(Enhanced chemiluminescent detection)
Western Blotting: 1:200-400(DAB detection)
Direct ELISA: 1:1000-5000(TMB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only