

FH Monoclonal Antibody

Catalog No. E-AB-22083

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:2L3
Purification	Protein A purification
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:1000-2000

Data



Western Blot analysis of 1) Mouse brain, 2) Rat brain using FH Monoclonal Antibody at dilution of 1:2000.

Observed Mw:50kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Fumarate hydratase,Also acts as a tumor suppressor

For Research Use Only