

EBFP Monoclonal Antibody

Catalog No. E-AB-22095

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:3C1
Purification	Protein A purification
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:3000-10000

Data



www.elabscience.com

Western Blot analysis of recombinant EBFP protein using EBFP Monoclonal Antibody at dilution of 1:5000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only