

## GPR17 Polyclonal Antibody

Catalog No. E-AB-18020

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                         |
|---------------------|---------------------------------------------------------|
| <b>Reactivity</b>   | Human, Mouse, Rat                                       |
| <b>Immunogen</b>    | Synthetic peptide of human GPR17                        |
| <b>Host</b>         | Rabbit                                                  |
| <b>Isotype</b>      | IgG                                                     |
| <b>Purification</b> | Antigen affinity purification                           |
| <b>Conjugation</b>  | Unconjugated                                            |
| <b>Buffer</b>       | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4 |

### Applications Recommended Dilution

|              |                |
|--------------|----------------|
| <b>IHC</b>   | 1:30-1:150     |
| <b>ELISA</b> | 1:5000-1:10000 |

### Data



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using GPR17 Polyclonal Antibody at dilution of 1:25(×200)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

GPR17 (G Protein-Coupled Receptor 17) is a Protein Coding gene. Among its related pathways are Signaling by GPCR and Nucleotide-like (purinergic) receptors. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and chemokine receptor activity. An important paralog of this gene is P2RY10.

### For Research Use Only