

GPR22 Polyclonal Antibody

Catalog No. E-AB-18022

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Synthetic peptide of human GPR22
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:40-1:200

Data




Western blot analysis of HepG2 cell lysate using GPR22 Polyclonal Antibody at dilution of 1:500

Observed Mw: Refer to figures
Calculated Mw: 49 kDa



Immunohistochemistry of paraffin-embedded Human cervical cancer tissue using GPR22 Polyclonal Antibody at dilution of 1:50 (x200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using GPR22 Polyclonal Antibody at dilution of 1:50 (x200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

This gene is a member of the G-protein coupled receptor 1 family and encodes a multi-pass membrane protein.