

## KBTBD11 Polyclonal Antibody

**Catalog No.** E-AB-18161

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Human, Mouse
<b>Immunogen</b>	Synthetic peptide of human KBTBD11
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol,pH7.4

### Applications Recommended Dilution

**WB** 1:1000-1:5000

### Data



Western blot analysis of Human metastatic papillary carcinoma(thyroid cancer) tissue lysate using KBTBD11 Polyclonal Antibody at dilution of 1:2300

**Observed Mw:Refer to figures**  
**Calculated Mw:66 kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

The function of this protein remains unknown.

### For Research Use Only