

ZNHIT1 Polyclonal Antibody

Catalog No. E-AB-18329

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:30-1:150

Data




Western blot analysis of 293T cell lysate using ZNHIT1 Polyclonal Antibody at dilution of 1:400

Observed Mw: Refer to figures
Calculated Mw: 18 kDa



Immunohistochemistry of paraffin-embedded Human tonsil tissue using ZNHIT1 Polyclonal Antibody at dilution of 1:40 (x200)



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ZNHIT1 Polyclonal Antibody at dilution of 1:40 (x200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

ZNHIT1, Zinc finger HIT domain containing protein 1, appears to play a role in p53-mediated apoptosis induction.

For Research Use Only