

NDUFB9 Polyclonal Antibody

Catalog No. E-AB-18339

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human NDUFB9
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of Hela cell and human fetal liver tissue A549 and Jurkat cell using NDUFB9 Polyclonal Antibody at dilution of 1:400

Observed Mw: Refer to figures
Calculated Mw: 22 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9 is an enzyme that in humans is encoded by the NDUFB9 gene. NADH dehydrogenase (ubiquinone) 1 beta subcomplex subunit 9 is an accessory subunit of the NADH dehydrogenase (ubiquinone) complex, located in the mitochondrial inner membrane. It is also known as Complex I and is the largest of the five complexes of the electron transport chain.

For Research Use Only