

S100A16 Polyclonal Antibody

Catalog No. E-AB-18345

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB 1:500-1:2000,
ELISA
1:5000-1:10000

Data



Western blot analysis of MCF-7 PC-3 and A431 cell
using S100A16 Polyclonal Antibody at dilution of
1:250

Observed Mw: Refer to figures
Calculated Mw: 12 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

S100A16 (S100 Calcium Binding Protein A16) is a Protein Coding gene. Diseases associated with S100A16 include Robinow Syndrome, Autosomal Recessive and Robinow Syndrome. Gene Ontology (GO) annotations related to this gene include protein homodimerization activity.

For Research Use Only