

MMGT1 Polyclonal Antibody

Catalog No. E-AB-18353

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB 1:500-1:2000,
ELISA
1:5000-1:10000

Data



Western blot analysis of 293T and Hela cell using MMGT1 Polyclonal Antibody at dilution of 1:250

Observed Mw: Refer to figures
Calculated Mw: 15 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

MMGT1 (membrane magnesium transporter 1), also known as TMEM32 (transmembrane protein 32), is a 131 amino acid single-pass type I membrane protein of the endoplasmic reticulum, golgi apparatus and early endosome that belongs to the membrane magnesium transporter family. As such, MMGT1 mediates magnesium transport and exists as two alternatively spliced isoforms. The gene encoding MMGT1 maps to human chromosome Xq26.3. The X and Y chromosomes are the human sex chromosomes. Chromosome X consists of about 153 million base pairs and nearly 1,000 genes. The combination of a X and Y chromosome lead to normal male development while two copies of X lead to normal female development.

For Research Use Only