

ZNF248 Polyclonal Antibody

Catalog No. E-AB-18384

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:20-1:100

Data




Western blot analysis of Human kidney tissue using ZNF248 Polyclonal Antibody at dilution of 1:200

Observed Mw: Refer to figures
Calculated Mw: 67 kDa



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ZNF248 Polyclonal Antibody at dilution of 1:25 (×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using ZNF248 Polyclonal Antibody at dilution of 1:25 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

ZNF248 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 8 C2H2-type zinc fingers and 1 KRAB domain. ZNF248 may be involved in transcriptional regulation.

For Research Use Only