

ABHD14B Polyclonal Antibody

Catalog No. E-AB-18410

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:25-1:100

Data




Western blot analysis of LO2 cell and Human testis tissue lysates using ABHD14B Polyclonal Antibody at dilution of 1:350

Observed Mw: Refer to figures
Calculated Mw: 22 kDa



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using ABHD14B Polyclonal Antibody at dilution of 1:30 (×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using ABHD14B Polyclonal Antibody at dilution of 1:30 (×200)

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

ABHD14B is a 210 amino acid protein that localizes to both the cytoplasm and the nucleus where it exists as two alternatively spliced isoforms. Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

For Research Use Only