

ZNF22 Polyclonal Antibody

Catalog No. E-AB-18423

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human ZNF22
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:20-1:100

Data



Western blot analysis of Jurkat cell using ZNF22 Polyclonal Antibody at dilution of 1:400
Observed Mw: Refer to figures
Calculated Mw: 26 kDa



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ZNF22 Polyclonal Antibody at dilution of 1:25 (x200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Binds DNA through the consensus sequence 5'-CAATG-3'. May be involved in transcriptional regulation and may play a role in tooth formation (By similarity).

For Research Use Only