

ZNF24 Polyclonal Antibody

Catalog No. E-AB-18455

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human ZNF24
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:25-1:100

Data




Western blot analysis of HUVEC and Hela cell lysates using ZNF24 Polyclonal Antibody at dilution of 1:350

Observed Mw: Refer to figures
Calculated Mw: 42 kDa



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using ZNF24 Polyclonal Antibody at dilution of 1:35 (×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ZNF24 Polyclonal Antibody at dilution of 1:35 (×200)

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

ZNF24 (Zinc Finger Protein 24) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include DNA binding transcription factor activity and sequence-specific DNA binding. An important paralog of this gene is ZNF396.

For Research Use Only