

UQCRQ Polyclonal Antibody

Catalog No. E-AB-18511

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Rat
Immunogen	Full length fusion protein
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:30-1:150

Data



Western blot analysis of 293T cell lysates using UQCRQ Polyclonal Antibody at dilution of 1:450

Observed Mw: Refer to figures
Calculated Mw: 10 kDa



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using UQCRQ Polyclonal Antibody at dilution of 1:45 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes a ubiquinone-binding protein of low molecular mass. This protein is a small core-associated protein and a subunit of ubiquinol-cytochrome c reductase complex III, which is part of the mitochondrial respiratory chain.

For Research Use Only