

ATL3 Polyclonal Antibody

Catalog No. E-AB-18527

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human ATL3
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:50-1:300

Data



Western blot analysis of EPG2 and Jurkat cell lysates using ATL3 Polyclonal Antibody at dilution of 1:800

Observed Mw: Refer to figures
Calculated Mw: 61 kDa



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ATL3 Polyclonal Antibody at dilution of 1:85 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes a member of a family of dynamin-like, integral membrane GTPases. The encoded protein is required for the proper formation of the network of interconnected tubules of the endoplasmic reticulum. Mutations in this gene may be associated with hereditary sensory neuropathy type IF. Alternatively spliced transcript variants that encode distinct isoforms have been described.

For Research Use Only