

GATD3B Polyclonal Antibody

Catalog No. E-AB-18569

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human GATD3B
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:30-1:150

Data



Western blot analysis of 293T A549 HeLa K562 and HEPG2 cell lysates using GATD3B Polyclonal Antibody at dilution of 1:400
Observed Mw:Refer to figures
Calculated Mw:28 kDa



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using GATD3B Polyclonal Antibody at dilution of 1:35(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

GATD3B (Glutamine Amidotransferase Like Class 1 Domain Containing 3B) is a Protein Coding gene. Among its related pathways are Gastric Cancer Network 1. An important paralog of this gene is GATD3A.

For Research Use Only