

MAGOHB Polyclonal Antibody

Catalog No. E-AB-18776

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human MAGOHB
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:50-1:300


Data



Western blot analysis of Human prostate tissue Jurkat cell Hela and HL-60 cell lysates using MAGOHB Polyclonal Antibody at dilution of 1:400
Observed Mw: Refer to figures
Calculated Mw: 17 kDa



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using MAGOHB Polyclonal Antibody at dilution of 1:70(×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using MAGOHB Polyclonal Antibody at dilution of 1:70(×200)

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Involved in mRNA splicing and in the nonsense-mediated decay (NMD) pathway.

For Research Use Only