

NIPSNAP3A Polyclonal Antibody

Catalog No. E-AB-18802

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human NIPSNAP3A
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of Mouse brain tissue lysate
using NIPSNAP3A Polyclonal Antibody at dilution of
1:550

Observed Mw: Refer to figures
Calculated Mw: 28 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport.

For Research Use Only