

AARSD1 Polyclonal Antibody

Catalog No. E-AB-18843

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human AARSD1
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB 1:500-1:2000,
ELISA
1:5000-1:10000

Data



Western blot analysis of HeLa and A549 cell lysates using AARSD1 Polyclonal Antibody at dilution of 1:350

Observed Mw: Refer to figures
Calculated Mw: 45 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

AARSD1 (Alanyl-TRNA Synthetase Domain Containing 1) is a Protein Coding gene. Diseases associated with AARSD1 include Ovary Adenocarcinoma. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and aminoacyl-tRNA editing activity. An important paralog of this gene is PTGES3L-AARSD1.

For Research Use Only