

Phospho-HSP27 (Ser82) Polyclonal Antibody

Catalog No. E-AB-20897

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from human HSP27 around the phosphorylation site of Ser82
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:100-1:300

Data



Western Blot analysis of HuvEc cells with Phospho-HSP27 (Ser82) Polyclonal Antibody at dilution of 1:1000

Observed Mw:27kDa
Calculated Mw:23kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Involved in stress resistance and actin organization.

For Research Use Only