

IPO11 Polyclonal Antibody

Catalog No. E-AB-18875

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human IPO11
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:50-1:300

Data




Western blot analysis of K562 cell Human fetal brain tissue lysates using IPO11 Polyclonal Antibody at dilution of 1:500

Observed Mw: Refer to figures
Calculated Mw: 113 kDa



Immunohistochemistry of paraffin-embedded Human tonsil tissue using IPO11 Polyclonal Antibody at dilution of 1:70 (x200)



Immunohistochemistry of paraffin-embedded Human breast cancer tissue using IPO11 Polyclonal Antibody at dilution of 1:70 (x200)

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).