

UBD Polyclonal Antibody

Catalog No. E-AB-18910

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human UBD
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

IHC	1:50-1:300
------------	------------

Data



Immunohistochemistry of paraffin-embedded Human tonsil tissue using UBD Polyclonal Antibody at dilution of 1:50(×200)



Immunohistochemistry of paraffin-embedded Human lung cancer tissue using UBD Polyclonal Antibody at dilution of 1:50(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

UBD (Ubiquitin D) is a Protein Coding gene. Diseases associated with UBD include Hepatocellular Carcinoma. Among its related pathways are TRAF Pathway and fMLP Pathway. GO annotations related to this gene include proteasome binding. An important paralog of this gene is UBB.

For Research Use Only