

## ZNF581 Polyclonal Antibody

**Catalog No.** E-AB-18932

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                         |
|---------------------|---------------------------------------------------------|
| <b>Reactivity</b>   | Human                                                   |
| <b>Immunogen</b>    | Fusion protein of human ZNF581                          |
| <b>Host</b>         | Rabbit                                                  |
| <b>Isotype</b>      | IgG                                                     |
| <b>Purification</b> | Antigen affinity purification                           |
| <b>Conjugation</b>  | Unconjugated                                            |
| <b>Buffer</b>       | PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol, pH7.4 |

### Applications Recommended Dilution

|            |            |
|------------|------------|
| <b>IHC</b> | 1:50-1:300 |
|------------|------------|

### Data



Immunohistochemistry of paraffin-embedded Human tonsil tissue using ZNF581 Polyclonal Antibody at dilution of 1:50(×200)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

ZNF581 (Zinc Finger Protein 581) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include nucleic acid binding. An important paralog of this gene is ZNF524.

### For Research Use Only