

REG1B Polyclonal Antibody

Catalog No. E-AB-18985

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human REG1B
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:1000-1:5000
IHC	1:50-1:300

Data



Western blot analysis of Human pancreas tissues lysate using REG1B Polyclonal Antibody at dilution of 1:800

Observed Mw: Refer to figures
Calculated Mw: 19 kDa



Immunohistochemistry of paraffin-embedded Human tonsil tissue using REG1B Polyclonal Antibody at dilution of 1:55 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Might act as an inhibitor of spontaneous calcium carbonate precipitation. May be associated with neuronal sprouting in brain, and with brain and pancreas regeneration.

For Research Use Only